sagemax[®]

Instructions for Use

Paint 3D

sagemax.com

[en]	Instructions for Use Stains and Glazes	[cs]	Návod k použití Barvy a glazury
[de]	Gebrauchsinformation Malfarben und Glasuren	[sk]	Návod na používanie Farbivá a glazúry
[fr]	Mode d'emploi Maquillants et glasures	[hu]	Használati utasítás Karakterizáló festék és Glazúr
[it]	Istruzioni d'uso Supercolori e glasure	[sr]	Упутство за употребу Боје и глазури
[es]	Instrucciones de uso Maquillajes y glaseados	[mk]	Упатство за употреба Бои и глазури
[pt]	Instruções de Uso Pigmentos e Glazes	[bg]	Инструкции за употреба Оцветители и глазури
[sv]	Bruksanvisning Stains och glaze	[sq]	Udhëzimet e përdorimit Ngjyrues dhe Glazura
[da]	Brugsanvisning Malefarver og glasurer	[ro]	Instrucțiuni de utilizare Machiaje și glazuri
[fi]	Käyttöohjeet Maalivärit ja lasitteet	[uk]	Інструкція із застосування фарби і глазур
[no]	Bruksanvisning Malefarger og glasurer	[et]	Kasutusjuhend Värvid ja glasuurid
[nl]	Gebruiksaanwijzing Stains en Glazes	[lv]	Lietošanas instrukcija Iekrāsojumi un glazūras
[el]	Οδηγίες Χρήσης Χρωστικές και εφυσάλωματα	[lt]	Naudojimo instrukcija Dažai ir glazūra
[tr]	Kullanım Talimatları Boya ve Glazürler	[ar]	تعليمات الاستخدام الصبغات ومواد الصقل
[ru]	Инструкция по применению Красители и глазури	[zh-tw]	使用說明 染色劑與釉料
[pl]	Instrukcja stosowania Farby i glazury	[zh-cn]	使用说明 染色剂与釉料
[sl]	Navodila za uporabo Barve in glazure	[ja]	使用説明書 ステインとグレーズ
[hr]	Upute za uporabu Boje i Glazure	[ko]	사용 설명서 착색 및 글레이즈

Intended use

Intended purpose

Characterization and glazing of ceramic restorations in the anterior and posterior region in the CTE range of $9.4\text{--}17.5 \times 10^{-6}/\text{K}$ (25–500 °C).

Patient target group

- › Patients with permanent teeth

Intended users

- › Dentists (clinical procedure)
- › Dental laboratory technicians (fabrication of restorations in the dental laboratory)

Special training

None

Use

For dental use only.

Description

Paint 3D are universal staining, glazing and structuring pastes for the characterization and glazing of ceramic materials in the CTE range of $9.4\text{--}17.5 \times 10^{-6}/\text{K}$ (25–500 °C).

Product name	Product description
Paint 3D Gingiva Glaze	<ul style="list-style-type: none"> › Glazing of ceramic anterior and posterior restorations › Diluting of the following pastes: Paint 3D Gingiva
Paint 3D Glaze Fluo	Glazing of full-contour ceramic restorations in the anterior and posterior region
Paint 3D Dentin Paint 3D Intense Effect Paint 3D Enamel, Effect	Staining and characterization of ceramic restorations in the anterior and posterior region
Paint 3D Gingiva	Staining and characterization of the gingival area
Paint 3D Texture	Design and individualization of the morphology and surface structure of ceramic restorations in the anterior and posterior region, as well as the gingival area
Paint 3D Universal Liquid	Diluting of the following pastes: Paint 3D Dentin, Paint 3D Intense Effect, Paint 3D Enamel, Effect, Paint 3D Gingiva, Paint 3D Gingiva Glaze, Paint 3D Glaze Fluo

Technical specifications

Characteristics	Unit	Specifications*	
Flexural strength	MPa	≥ 50	
Chemical solubility	$\mu\text{g}/\text{cm}^2$	< 100	
Coefficient of thermal expansion (25–400°C)	$10^{-6}/\text{K}^{-1}$	Paint 3D Gingiva Glaze	9.3 ± 0.5
		Paint 3D Glaze Fluo	
		Paint 3D Dentin	9.1 ± 0.5
		Paint 3D Intense Effect	
		Paint 3D Enamel, Effect	
		Paint 3D Gingiva	
		Paint 3D Texture	9.2 ± 0.5
Glass transition temperature T_g	°C	460 ± 20	
Radioactivity U^{238}	Bq/g	≤ 1.0	

* The performance criteria defined by ISO 6872:2015 have been met. Classification: Dental ceramic type 1 / Class 1

Firing table

These firing parameters are guidance values.

Deviations (approx. $\pm 10\text{ }^{\circ}\text{C}/18\text{ }^{\circ}\text{F}$) may occur:

- Depending on the furnace generation
- Regional differences in the supply voltage
- Operating several electrical devices in the same electric circuit

Firing parameters

Stain and glaze firing when using the staining technique (monolithic) Paint 3D Glaze Fluo, Dentin, Enamel, Effect, Intense Effect, Glaze Gingiva, Gingiva Base, Gingiva Effect

Paint 3D	Stand-by temperature	Closing time	Heating rate	Firing temperature	Holding time	Vacuum on	Vacuum off	Long-term cooling
	B [$^{\circ}\text{C}/^{\circ}\text{F}$]	S* [min]	t [$^{\circ}\text{C}/^{\circ}\text{F}/\text{min}$]	T [$^{\circ}\text{C}/^{\circ}\text{F}$]	H [min]	V1 [$^{\circ}\text{C}/^{\circ}\text{F}$]	V2 [$^{\circ}\text{C}/^{\circ}\text{F}$]	L [$^{\circ}\text{C}/^{\circ}\text{F}$]
Zirconium oxide (ZrO_2) 1-4 units	403/757	IRT/6	45/81	710/1310	2:00	-	-	-
Zirconium oxide (ZrO_2) 5-7 units	403/757	IRT/6	40/72	720/1328	2:00	-	-	**
Zirconium oxide (ZrO_2) 8-14 units	403/757	IRT/6	30/54	730/1346	2:00	-	-	**
Lithium disilicate (LS_2) CAD	403/757	IRT/6	45/81	710/1310	2:00	-	-	***
Lithium disilicate (LS_2) Press	403/757	IRT/6	45/81	710/1310	2:00	-	-	-

*IRT standard mode

** In case of significant differences in the cross-sections of individual units within a restoration, long-term cooling L down to $500\text{ }^{\circ}\text{C}/932\text{ }^{\circ}\text{F}$ is recommended for each firing cycle.

*** For layer thicknesses exceeding 2 mm, long-term cooling L down to $500\text{ }^{\circ}\text{C}/932\text{ }^{\circ}\text{F}$ is required.

Texture firing in the texturing technique (monolithic) Paint 3D Texture / Paint 3D Texture Gingiva

Paint 3D	Stand-by temperature	Closing time	Heating rate	Firing temperature	Holding time	Vacuum on	Vacuum off	Long-term cooling
	B [$^{\circ}\text{C}/^{\circ}\text{F}$]	S* [min]	t [$^{\circ}\text{C}/^{\circ}\text{F}/\text{min}$]	T [$^{\circ}\text{C}/^{\circ}\text{F}$]	H [min]	V1 [$^{\circ}\text{C}/^{\circ}\text{F}$]	V2 [$^{\circ}\text{C}/^{\circ}\text{F}$]	L [$^{\circ}\text{C}/^{\circ}\text{F}$]
Zirconium oxide (ZrO_2) 1 unit	403/757	IRT/6	45/81	710/1310	2:00	550/1022	710/1310	-
Zirconium oxide (ZrO_2) 2-4 units	403/757	IRT/6	40/72	710/1310	2:00	550/1022	710/1310	-
Zirconium oxide (ZrO_2) 5-7 units	403/757	IRT/6	45/81	720/1328	2:00	550/1022	720/1328	500/932
Zirconium oxide (ZrO_2) 8-14 units	403/757	IRT/6	30/54	720/1328	2:00	550/1022	720/1328	500/932
Lithium disilicate (LS_2) CAD	403/757	IRT/6	45/81	710/1310	2:00	550/1022	710/1310	***
Lithium disilicate (LS_2) Press	403/757	IRT/6	45/81	710/1310	2:00	550/1022	710/1310	-

*IRT standard mode

*** For layer thicknesses exceeding 2 mm, long-term cooling L down to $500\text{ }^{\circ}\text{C}/932\text{ }^{\circ}\text{F}$ is required.

Stain and Glaze firing in the layering technique (full veneering / partial veneering) Paint 3D Glaze Fluo, Dentin, Enamel, Effect, Intense Effect, Glaze Gingiva, Gingiva Base, Gingiva Effect

Paint 3D	Stand-by temperature	Closing time	Heating rate	Firing temperature	Holding time	Vacuum on	Vacuum off	Long-term cooling
	B [°C/°F]	S* [min]	t [°C/°F/min]	T [°C/°F]	H [min]	V1 [°C/°F]	V2 [°C/°F]	L [°C/°F]
Veneering ceramics (all-ceramic) / zirconium oxide (ZrO ₂) 1-4 units	403/757	IRT/6	45/81	710/1310	1:00	-	-	450/842
Veneering ceramics (all-ceramic) / zirconium oxide (ZrO ₂) 5-7 units	403/757	IRT/6	40/72	720/1328	1:00	-	-	450/842
Veneering ceramics (all-ceramic) / zirconium oxide (ZrO ₂) 8-14 units	403/757	IRT/6	30/54	720/1328	1:00	-	-	450/842
Veneering ceramics (all-ceramic) / lithium disilicate (LS ₂) CAD	403/757	IRT/6	60/108	710/1310	1:00	-	-	-
Veneering ceramics (all-ceramic) / lithium disilicate (LS ₂) Press	403/757	IRT/6	60/108	710/1310	1:00	-	-	-
Veneering ceramics (metal-ceramic)	403/757	IRT/6	60/108	750/1382	1:00	-	-	-

*IRT standard mode

Indications

None

Areas of application:

Characterization and glazing of ceramic materials:

- › Zirconium oxide
- › Milling ceramics
- › Press ceramics
- › Veneering ceramics (all-ceramics & metal-ceramics)

Contraindications

- › The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Restrictions of use

The product must not be used in the following cases:

- › Staining and glazing of partially crystallized lithium disilicate restorations (in their blue state)
- › Reuse of the final restoration

Side effects

None known to date

Interactions

None known to date

Clinical benefit

- › Restoration of esthetics

Composition

- › Paint 3D Gingiva Glaze
Components: Alkali aluminosilicate glass, 1,3-butanediol, glycerine
- › Paint 3D Glaze Fluo
Components: Alkali aluminosilicate glass, 1,3-butanediol, glycerine
- › Paint 3D Dentin
Components: Alkali aluminosilicate glass, pigments, 1,3-butanediol, glycerine, water
- › Paint 3D Intense Effect
Components: Alkali aluminosilicate glass, pigments, 1,3-butanediol, glycerine, water
- › Paint 3D Enamel, Effect
Components: Alkali aluminosilicate glass, pigments, 1,3-butanediol, glycerine, water
- › Paint 3D Gingiva
Components: Alkali aluminosilicate glass, pigments, 1,3-butanediol, glycerine, water
- › Paint 3D Texture
Components: Alkali aluminosilicate glass, pigments, 1,3-butanediol, water
- › Paint 3D Universal Liquid
Components: 1,3-butanediol

Application

Preparing the ceramic surface

NOTICE! Depending on the material, the restoration surface must be prepared in varying ways. Observe the Instructions for Use of the respective material.

NOTICE! Smooth or prepolish the areas for which a higher gloss is desired after glaze firing using silicone polishers.

Monolithic restorations made of zirconium oxide (ZrO₂)

NOTICE! Make changes to sintered restorations only if absolutely necessary.

1. If absolutely necessary: Prepare the restoration with suitable grinding instruments (see Instructions for Use) at low pressure and low speed.
2. Optional: Blast the restoration with aluminium oxide (Al₂O₃) 25-70 µm at 1 bar/15 psi pressure or 70-110 µm at 1.5 bar/22 psi pressure.
3. Smooth the incisal and occlusal contact points and the basal side of the bridge connectors using rubber polishers.
4. Thoroughly clean the restoration with running water or the steam jet.
5. Dry the restoration.

Monolithic restorations made of lithium disilicate (LS₂) CAD

NOTICE! Do not blast lithium disilicate restorations with Al₂O₃ or polishing beads.

- ✓ The restoration is in its precrystallized (blue) state.
1. Prepare restoration using diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
 2. Clean the restoration with ultrasound in a water bath and / or with a steam jet.
 3. Dry the restoration.
 4. Crystallize the restoration.

Monolithic restorations made of lithium disilicate (LS₂) Press

1. Prepare the restoration using ceramic- and / or diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
2. Blast restoration with Al₂O₃, 100 µm at 1 bar/15 psi pressure.
3. Thoroughly clean the restoration with running water or the steam jet.
4. Dry the restoration.

Monolithic restorations made of leucite-reinforced glass-ceramic CAD

1. Prepare the restoration using ceramic- and / or diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
2. Carefully prepare the restoration margins using rubber and silicone polishers. Observe the Instructions for Use.
3. Thoroughly clean the restoration with running water or the steam jet.
4. Dry the restoration.

Veneering ceramics

1. Prepare the restoration using ceramic- and / or diamond-bonded grinding instruments at low pressure and low speed. Observe the Instructions for Use.
2. Thoroughly clean the restoration with running water or the steam jet.
3. Dry the restoration.

Characterizing and glazing the restoration

- ✓ The restoration is free of dirt and grease residues.
- ✓ The restoration demonstrates a smooth surface without sharp edges / ridges.
- ✓ Paint 3D was thoroughly mixed with a metal-free spatula.

NOTICE! Avoid pooling and too thick stain layers. More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers.

White / red esthetics for monolithic restorations

Wet-in-wet staining / one-shot

1. Apply undiluted Paint 3D Glaze Fluo on the teeth in an even layer using a brush.
2. Stain the teeth with Paint 3D Dentin/Enamel and adjust shade in accordance with the shade guide.
3. Characterize the teeth with Paint 3D Enamel, Effect and Paint 3D Intense Effect.
4. Optional: Apply undiluted Paint 3D Gingiva Base in the gingival area (self glazing).
5. Optional: Characterize the gingival area with Paint 3D Gingiva (self glazing).
6. Conduct Stain and Glaze firing. See firing parameters.
7. Repeat Steps 1-6, if required.

Conventional staining technique

1. Apply Paint 3D Glaze Fluo alone or slightly diluted with Paint 3D Universal Liquid on the teeth in a thin layer using a brush.
2. Stain the teeth with Paint 3D Dentin/Enamel and adjust shade in accordance with the shade guide.
3. Characterize the teeth with Paint 3D Enamel, Effect and Paint 3D Intense Effect.
4. Optional: Apply undiluted Paint 3D Gingiva Base in the gingival area (self glazing).
5. Conduct Stain firing. See firing parameters.
6. Repeat Steps 1-5, if required.
7. Apply undiluted Paint 3D Glaze Fluo on the teeth in an even layer using a brush.
8. Optional: Characterize the gingival area with Paint 3D Gingiva (self glazing).
9. Conduct Stain and Glaze firing. See firing parameters.
10. Repeat Steps 7-9, if required.

Texturing technique

NOTICE! For better sculptability of Paint 3D Texture moisten the brush / instrument with Paint 3D Universal Liquid. Paint 3D Texture must not be diluted with Paint 3D Universal Liquid .

NOTICE! After application of Paint 3D Texture check the contact area to the antagonist, as well as the contact points.

- ✓ Application in the wet-in-wet technique or the conventional technique has been performed.
- 1. Apply undiluted Paint 3D Texture to the glazed surface of the restoration (self glazing).
- 2. Sculpt the morphology of the restoration.
- 3. Incorporate the microtexture in the surface using instruments and brushes.
- 4. Conduct texture firing. See firing parameters.

White / red esthetics for veneering materials

Layering technique / one-shot

1. Apply Paint 3D Gingiva Glaze in a thin even layer to the teeth using a brush.
2. Stain the teeth with Paint 3D Dentin/Enamel and adjust shade in accordance with the shade guide.
3. Characterize the teeth with Paint 3D Enamel, Effect and Paint 3D Intense Effect.
4. Optional: Apply undiluted Paint 3D Gingiva Base in the gingival area (self glazing).
5. Optional: Characterize the gingival area with Paint 3D Gingiva and Paint 3D Intense Effect.
6. Conduct Stain and Glaze firing. See firing parameters.
7. Repeat Steps 1–6, if required.

Cementation

NOTICE! All working steps in the oral cavity of the patient must exclusively be performed by a qualified professional (e.g. dentist).

- ✓ The cementation surface of the restoration is free of ceramic material.
- 1. Cement the restoration. Observe the Instructions for Use.

Safety information

- › The product has been developed solely for use in dentistry. Process according to the Instructions for Use.
- › In the case of serious incidents related to the product, please contact DSSM AG, Im alten Riet 9, 9494 Schaan/Liechtenstein, vigilance@dssm.group and your responsible competent authority.
- › Do not inhale ceramic dust during finishing. Use suction equipment, a face mask and protective goggles.
- › If applied frequently or for a prolonged period of time, peroxide-based professional bleaching agents (carbamide peroxide; hydrogen peroxide) as well as acidulated phosphate fluorides used for caries prevention can make the surface of restorations rough and matte.
- › Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact DSSM AG or your local sales partner.

Supporting documents

Document	To be found:
Current version of Instructions for Use	dssm-eifu.sagemax.com
Structure of Instructions for Use and Warnings	dssm-eifu.sagemax.com
Explanation of symbols	dssm-eifu.sagemax.com
Safety Data Sheet (SDS)	www.sagemax.com
Summary of Safety and Clinical Performance - SSCP Basic UDI-DI: 763071662APA3D001DB	https://ec.europa.eu/tools/eudamed

Disposal information

Dispose contaminated products, remaining stocks or removed restorations according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- › Chipping / fracture of restorative material may lead to swallowing or inhalation of material and to repeated dental treatment.

Shelf life and storage

- › Storage temperature 2-28 °C (36-82 °F)
- › Store the product free of vibration.
- › Store the product in a dry place.
- › Protect the product from sunlight.
- › Do not use the product after the indicated expiry date.
- › Date of expiration: See note on packaging

Additional information

Keep material out of the reach of children!

The manufacturer accepts no liability for damages resulting from improper use of the product. Furthermore, the user is responsible for determining the product's suitability for the user's own purposes before use.

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